

ARAGEL®

Ophthalmic gel based on Hyaluronic Acid sodium salt, Arabinogalactan and Trehalose
Preservative-free

DESCRIPTION

Lubricant ophthalmic gel, based on Hyaluronic Acid sodium salt, Arabinogalactan and Trehalose, preservative free, it is suitable for both day and night use.

The Hyaluronic Acid sodium salt gives a high viscoelasticity to the product, ensuring a more effective hydration to the ocular surface, such that dry eye symptoms are alleviated.

The Arabinogalactan is a natural mucopolysaccharide extracted from the bark of the Canadian Larch tree, it has lubricating properties, and it is effective in corneal protection by restoring the mucin layer of the corneal epithelium.

The Trehalose, in combination with Vitamin E TPGS, gives the product antioxidant properties, as well as high moisturizing properties.

The combination of Arabinogalactan, Trehalose and Hyaluronic Acid sodium salt confers to ARAGEL® a long-lasting protective and moisturizing activity that, in combination with the emollient action of D-Panthenol, produces rapid relief.

The addition of Poloxamer 407 gives consistency of Gel in contact with the ocular surface due to a higher temperature.

INDICATIONS

ARAGEL® is a lubricating ophthalmic gel indicated in the treatment of moderate-severe dry eye syndrome, for foreign body sensation and for protection, hydration and lubrication of the ocular surface after corneal surface surgery.





The product can also be used during contact lens wear.

FREQUENCY OF USE

Instill 1 drop of ARAGEL® in each eye, 2 times a day or as advised by your ophthalmologist.

INSTRUCTIONS FOR USE

1. Open the protective foil pouch
2. Carefully clean hands before the use
3. Separate the single-dose vial from the strip and open it by twisting the cap
4. Gently press the single-dose vial and instill 1 drop per eye
5. After the instillation, close the vial and reuse the residue solution within 12 hours
6. Store the unused vials in the original foil pouch
7. The product is single patient to avoid possible contamination
8. The product can be used with contact lenses



COMPOSITION

Arabinogalactan, Trehalose, Hyaluronic acid sodium salt, Vitamin E-TPGS, D-Panthenol, Poloxamer 407, boric acid, sodium tetraborate decahydrate, purified water.

PACKAGING

The box is composed of 2 foil pouches each one containing one strip of 5 single-dose vials. Every single-dose vial is filled at 0,5 ml and it is reclosable.

The box is composed of 3 foil pouches each one containing one strip of 5 single-dose vials. Every single-dose vial is filled at 0,5 ml and it is reclosable.

The box is composed of 4 foil pouches each one containing one strip of 5 single-dose vials. Every single-dose vial is filled at 0,5 ml and it is reclosable.

The box is composed of 5 foil pouches each one containing one strip of 5 single-dose vials. Every single-dose vial is filled at 0,5 ml and it is reclosable.

The box is composed of 6 foil pouches each one containing one strip of 5 single-dose vials. Every single-dose vial is filled at 0,5 ml and it is reclosable.

WARNINGS

- The product is for ophthalmic use only, do not ingest or use otherwise than the indication of use
- Once open, the single-dose vial can be reclosed with its cap and the residue solution can be reused within 12 hours to avoid possible contaminations of the product
- Keep out of reach of children
- During the instillation do not touch the eye or other surfaces with the nozzle of the single-dose vial
- Do not use the product in case of known incompatibility to even one single component of the formulation
- If necessary, check with your optician for the compatibility of the ophthalmic gel with the contact lenses use
- In case of adverse reaction during the use, suspend the treatment and consult your doctor or pharmacist.
- Consult your doctor for use during pregnancy, lactation and for children
- The treatment with ophthalmic drugs or other ophthalmic medical devices must be carried out under the supervision of a doctor who can evaluate possible interactions.
- In case of compromised ocular surface or damage cornea consult your doctor
- Do not use the product after the expiration date indicated on the package
- In the event of any serious claims or incidents concerning the device, it should be reported the incident to the manufacturer and the competent authority of the Member State in which the user is established, to provide more information on the safety of this medical device.

Subject:

ARAGEL® – Single-dose ophthalmic gel


Rev.01 of January 2025

Pages 3/3

VALIDITY: 36 months with the package intact



: Enable Innovations S.r.l. - Via Claudio Tolomeo n°26; 41012 Carpi (MO), Italy



Rev.01 January 2025

Optional information:



: XXXXXXX